

Bathboard

A guide for total hip joint replacement patients and their carers

Your Occupational Therapist is

 9288 2211 or pager _____

Setting up the bathboard

Your bathboard will need to be installed by a family member/friend or an equipment company. Your Occupational Therapist (OT) will show your carer how to do this if required.

1. Make sure there is a lip of 4-5cm or more on both edges of the bath.
2. Place the bathboard across the bath so the board sits on both edges. Tighten the clamps underneath to secure.
3. Where possible, the bathboard should be placed under the stream of water and where it is easy to reach the taps.
4. Check with your OT that the height of the bathboard is appropriate for meeting your hip precautions.



Getting in the bath

1. Ensure the bathboard is secure.
2. Back up to the bathboard until you can feel the back of your legs against the edge of the bath.
3. Reach back for the board with both hands.

4. Slide your operated leg forward and sit down on the bathboard.



5. Shuffle/move your bottom back on the bathboard while leaning back.



6. Turn (without twisting) and lift both legs into the bath while leaning back slightly.



Bathboard

A guide for total hip joint replacement patients and their carers

Getting out of the bath

1. Lean back slightly and turn (without twisting) on the bathboard and lift both legs outside the bath.



2. Slide your operated leg out in front. Push up from bathboard using both hands. Do not bend too much at the hip.



If you have any questions regarding this equipment or task, please contact your Occupational Therapist.